ARJOHUNTLEIGH GETINGE GROUP





This fixed-height keyhole shaped tub provides an ideal solution when a care facility requires an advanced sit bath designed for good immersion.

The sit bath is a space-saving alternative that is regarded as a traditional bathing method in some parts of the world. *Century*[™] makes the ArjoHuntleigh bathing system range fully comprehensive by offering the best available sit bath. *Century* enables smoother transfers and it's a keyhole-shape gives the carer easier access to the patient during bathing sessions.



...with people in mind


IMPROVED SIT BATH SYSTEM



Century integrates with *Alenti* and *Calypso.*



Two sensors scalding protection and Auto-fill, one level with possible afterfilling. The electronic control functions are operated by touch pad



Digital display of temperature for bath water and fill/shower.




Built-in disinfection shower handle and automatic disinfection of the hydromassage system.





Hydromassage generates a gentle massaging effect that stimulates peripheral blood circulation and eases muscle tension.



Hydrosound offers a gentle touch for cleaning sensitive skin. An integrated ultrasound unit in the tub generates millions of jet streams.



Product information	
Maximum overall width	850 mm (33 ½")
Maximum overall length	1610 mm (63 ¼")
Maximum tub edge height	860 mm (33 ¾")
Total weight filled with water	420 kg/925 lbs
Tub volume	200 L (53 gal US)
For more information of Contury, please visit	

For more information of Century, please visit www.ArjoHuntleigh.com

ArjoHuntleigh is a branch of Arjo Ltd Med. AB. Only ArjoHuntleigh designed parts, which are designed specifically for the purpose, should be used on the equipment and products supplied by ArjoHuntleigh. As our policy is one of continuous development we reserve the right to modify designs and specifications without prior notice. @ and ™ are trademarks belonging to the ArjoHuntleigh group of companies. © ArjoHuntleigh, 2010

ArjoHuntleigh

310-312, Dallow Road Luton Bedfordshire LU1 1TD United Kingdom • Phone: +44 1582 745 800 • Fax: +44 1582 745 866 E-mail: international@arjohuntleigh.com • www.ArjoHuntleigh.com